



**By AI**

# Intracranial Hemorrhage Detection

**Yuming Chen**

**Bo Feng**



**兰州大学**  
LANZHOU UNIVERSITY




# Project Overview

# Overview

- **Serious** Health Problem
- **Rapid & Intensive** Treatment
- **Complicated** Process




- **Detection Subtypes(5)**
- **Multi**-class Classification
- InceptionNetV3 & Spark



# Data Information

# Data Origin



Featured Prediction Competition

## RSNA Intracranial Hemorrhage Detection

Identify acute intracranial hemorrhage and its subtypes

**\$25,000**  
Prize Money

RSNA Radiological Society of North America · 1,345 teams · 2 years ago

Overview Data Code Discussion Leaderboard Rules Team My Submissions **Late Submission**

Overview

**Description**

**Hemorrhage Types**

**Evaluation**

Intracranial hemorrhage, bleeding that occurs inside the cranium, is a serious health problem requiring rapid and often intensive medical treatment. For example, intracranial hemorrhages account for approximately 10% of strokes in the U.S., where stroke is the fifth-leading cause of death. Identifying the location and type of any hemorrhage present is a critical step in treating the patient.

and often time consuming.

In this competition, your challenge is to build an algorithm to detect acute intracranial hemorrhage and its subtypes.

You'll develop your solution using a rich image dataset provided by the Radiological Society of North America (RSNA®) in collaboration with members of the American Society of Neuroradiology and MD.ai.

If successful, you'll help the medical community identify the presence, location and type of hemorrhage in order to quickly and effectively treat affected patients.

Challenge participants may be invited to present their AI models and methodologies during an award ceremony at the RSNA Annual Meeting which will be held in Chicago, Illinois, USA, from December 1-6, 2019.


**Collaborators**

Four research institutions provided large volumes of de-identified CT studies that were assembled to create the challenge dataset: Stanford University, Thomas Jefferson University, Unity Health Toronto and Universidade Federal de São Paulo (UNIFESP). The American Society of Neuroradiology (ASNR) organized a cadre of more than 60 volunteers to label over 25,000 exams for the challenge dataset. ASNR is the world's leading organization for the future of neuroradiology representing more than 5,300 radiologists, researchers, interventionalists, and imaging scientists. MD.ai provided tooling and support for the data annotation process.

The RSNA is an international society of radiologists, medical physicists and other medical professionals with more than 54,000 members from 146 countries across the globe. They see the potential for AI to assist in detection and classification of hemorrhages in order to prioritize and expedite their clinical work.

A full set of acknowledgments can be found on this page.

kaggle [www.kaggle.com](https://www.kaggle.com)



Search News About Donate Login

RSNA®  
Radiological Society  
of North America

Membership Annual Meeting Journals Education Research Practice Tools **COVID-19**

**RSNA 2021**  
REDEFINING RADIOLOGY  
November 28 to December 2

**Join us in Chicago for RSNA 2021!**

Reserve your hotel now for the best pricing, selection and flexibility.  
Meeting registration opens July 21.

[Reserve now](#)

**INNOVATION**

**We empower radiology's innovators**

For more than 100 years we've been at the forefront of our field, convening thought leaders at our annual meeting and developing tools to elevate your practice.

RSNA® [www.rsna.org](https://www.rsna.org)




# Data Format

## CSV

id	label
ID_52c205043	epidural
ID_15e6dc40a	intraventricular
ID_01a97add4	intraparenchymal
ID_9fed7e3db	subdural
ID_903d2ef48	subdural
ID_e0942808d	epidural
ID_d4d31f695	subdural
.....	
ID_c64f8db24	subarachnoid
ID_d202babf2	subarachnoid
ID_9ffd25faa	subdural

19877

## DICOM



0.53MB

512X512

```
Dataset.file_meta
(0002, 0000) File Meta Information Group Length  UL: 176
(0002, 0001) File Meta Information Version       OB: b'\x00\x01'
(0002, 0002) Media Storage SOP Class UID        UI: CT Image Storage
(0002, 0003) Media Storage SOP Instance UID     UI: 9999.1345638575715201386455629118010580
23921
(0002, 0010) Transfer Syntax UID                UI: Explicit VR Little Endian
(0002, 0012) Implementation Class UID          UI: 1.2.40.0.13.1.1.1
(0002, 0013) Implementation Version Name       SH: 'dcm4che-1.4.38'

(0008, 0018) SOP Instance UID                  UI: ID_7ce5ae372
(0008, 0060) Modality                          CS: 'CT'
(0010, 0020) Patient ID                        LO: 'ID_314d3781'
(0020, 000d) Study Instance UID                UI: ID_af0d670f5a
(0020, 000e) Series Instance UID              UI: ID_4195a174b7
(0020, 0010) Study ID                          SH: ''
(0020, 0032) Image Position (Patient)          DS: [-126.408875, -126.408875, -244.165497]
(0020, 0037) Image Orientation (Patient)       DS: [1.000000, 0.000000, 0.000000, 0.000000
0, 1.000000, 0.000000]
(0028, 0002) Samples per Pixel                 US: 1
(0028, 0004) Photometric Interpretation        CS: 'MONOCHROME2'
(0028, 0010) Rows                             US: 512
(0028, 0011) Columns                          US: 512
(0028, 0030) Pixel Spacing                     DS: [0.494750976563, 0.494750976563]
(0028, 0100) Bits Allocated                   US: 16
(0028, 0101) Bits Stored                      US: 16
(0028, 0102) High Bit                         US: 15
(0028, 0103) Pixel Representation              US: 1
(0028, 1050) Window Center                     DS: "35.0"
(0028, 1051) Window Width                     DS: "135.0"
(0028, 1052) Rescale Intercept                 DS: "-1024.0"
(0028, 1053) Rescale Slope                    DS: "1.0"
(7fe0, 0010) Pixel Data                        OW: Array of 524288 elements
```

10434.73MB



## Pydicom

[gitter](#) [join chat](#) [🔗](#)

Dicom (Digital Imaging in Medicine) is the bread and butter of medical image datasets, storage and transfer. This is the future home of the Pydicom documentation. If you are a Python developer looking to get started with Dicom and Python, this will be the place to learn and contribute! For now, here are some helpful links, and general plan for some of the code bases in the organization. If you want to come and chat, find our community on Gitter, or post an issue on one of our repos.

### Modules

#### Pydicom

If you want to work with dicom datasets, you should use [pydicom](#) [🔗](#). We have started a base of docs here, and see the [documentation](#) [🔗](#) for you to get started.

#### Pydicom

Information>

Software >

Datasets >

Containers>



[pydicom.github.io](https://pydicom.github.io)




The screenshot shows the DICOM website homepage. At the top, there is a navigation bar with the DICOM logo, a search bar, and a login link. Below the navigation bar is a large hero section with a blue background featuring a globe and network lines. The hero section contains the text "DICOM makes medical imaging information interoperable" and a "Learn More" button. Below the hero section are three smaller sections: "DICOM Standard", "Working Groups", and "Presentations", each with a representative image.



[www.dicomstandard.org](http://www.dicomstandard.org)


	Intraparenchymal	Intraventricular	Subarachnoid	Subdural	Epidural
<b>Location</b>	Inside of the brain	Inside of the ventricle	Between the arachnoid and the pia mater	Between the Dura and the arachnoid	Between the dura and the skull
<b>Imaging</b>					
<b>Mechanism</b>	High blood pressure, trauma, arteriovenous malformation, tumor, etc	Can be associated with both intraparenchymal and subarachnoid hemorrhages	Rupture of aneurysms or arteriovenous malformations or trauma	Trauma	Trauma or after surgery
<b>Source</b>	Arterial or venous	Arterial or venous	Predominantly arterial	Venous (bridging veins)	Arterial
<b>Shape</b>	Typically rounded	Conforms to ventricular shape	Tracks along the sulci and fissures	Crescent	Lentiform
<b>Presentation</b>	Acute (sudden onset of headache, nausea, vomiting)	Acute (sudden onset of headache, nausea, vomiting)	Acute (worst headache of life)	May be insidious (worsening headache)	Acute (skull fracture and altered mental status)






# **Exploratory Data Analysis**

label	count
subarachnoid	4177
epidural	3145
subdural	4132
intraventricular	4072
intraparenchymal	4351




# Visualization for Different Windows




**Window:** also known as 'grey-level mapping', contrast stretching, histogram modification or contrast enhancement' is the process in which the CT image greyscale component of an image is manipulated via the CT numbers; doing this will change the appearance of the picture to highlight particular structures. The brightness of the image is adjusted via the window level. The contrast is adjusted via the window width.

# Visualization for Different Windows


No Windows




Brain Windows



Three Channels Windows




Brain + Subdural + Bone Windows






# Visualization for Different Classes


Instances with Intraparenchymal Hemorrhage




Instances with Epidural Hemorrhage




Instances with Intraventricular Hemorrhage



Subarachnoid




Subdural



# Data Preprocess

# Data Preprocess

## Get windowed image



Rescale image into **299X299**

Apply operation by RDD method

## Add Label

id	image	label
ID_7ce5ae372	0.0, 0.0, 0.0, 0.0, 0.0...	intraparenchymal
ID_ace564fe2	0.0, 0.0, 0.0, 0.0, 0.0...	intraventricular
ID_9416259b1	0.0, 0.0, 0.0, 0.0, 0.0...	epidural
ID_928806565	0.0, 0.0, 0.0, 0.0, 0.0...	subarachnoid
ID_e97dd15eb	0.0, 0.0, 0.0, 0.0, 0.0...	intraparenchymal

only showing top 5 rows

## Flatten

id	image
ID_7ce5ae372	0.0, 0.0, 0.0, 0.0, 0.0...
ID_ace564fe2	0.0, 0.0, 0.0, 0.0, 0.0...
ID_9416259b1	0.0, 0.0, 0.0, 0.0, 0.0...
ID_928806565	0.0, 0.0, 0.0, 0.0, 0.0...
ID_e97dd15eb	0.0, 0.0, 0.0, 0.0, 0.0...

only showing top 5 rows

## Map Label in to Numeric Value

label	labelIndex
subarachnoid	3.0
intraventricular	2.0
epidural	0.0
intraparenchymal	1.0
subdural	4.0


## Split

Size of training dataset: 15951  
Size of testing dataset: 3926

Save

**Model**






[pytorch.org](https://pytorch.org)

## SparkTorch

build passing license MIT pypi v0.1.2



This is an implementation of Pytorch on Apache Spark. The goal of this library is to provide a simple, understandable interface in distributing the training of your Pytorch model on Spark. With SparkTorch, you can easily integrate your deep learning model with a ML Spark Pipeline. Underneath the hood, SparkTorch offers two distributed training approaches through tree reductions and a parameter server. Through the api, the user can specify the style of training, whether that is distributed synchronous or hogwild.

### Why should I use this?



Like SparkFlow, SparkTorch's main objective is to seamlessly work with Spark's ML Pipelines. This library provides three core components:

- Data parallel distributed training for large datasets. SparkTorch offers distributed synchronous and asynchronous training methodologies. This is useful for training very large datasets that do not fit into a single machine.
- Full integration with Spark's ML library. This ensures that you can save and load pipelines with your trained model.
- Inference. With SparkTorch, you can load your existing trained model and run inference on billions of records in parallel.

On top of these features, SparkTorch can utilize barrier execution, ensuring that all executors run concurrently during training (This is required for synchronous training approaches).



[github.com/dmmiller612/sparktorch](https://github.com/dmmiller612/sparktorch)



Model Structure

```
class InceptionNet3(nn.Module):
    def __init__(self, num_classes=5):
        super(InceptionNet3, self).__init__()
        self.Conv2d_1a_3x3 = BasicConv2d(3, 32, kernel_size=
3, stride=2)
        self.Conv2d_2a_3x3 = BasicConv2d(32, 32, kernel_size
=3)
        self.Conv2d_2b_3x3 = BasicConv2d(32, 64, kernel_size
=3, padding=1)
        self.Conv2d_3b_1x1 = BasicConv2d(64, 80, kernel_size
=1)
        self.Conv2d_4a_3x3 = BasicConv2d(80, 192, kernel_siz
e=3)
        self.Mixed_5b = InceptionA(192, pool_features=32)
        self.Mixed_5c = InceptionA(256, pool_features=64)
        self.Mixed_5d = InceptionA(288, pool_features=64)
        self.Mixed_6a = InceptionB(288)
        self.Mixed_6b = InceptionC(768, channels_7x7=128)
        self.Mixed_6c = InceptionC(768, channels_7x7=160)
        self.Mixed_6d = InceptionC(768, channels_7x7=160)
        self.Mixed_6e = InceptionC(768, channels_7x7=192)
        self.Mixed_7a = InceptionD(768)
        self.Mixed_7b = InceptionE(1280)
        self.Mixed_7c = InceptionE(2048)
        self.fc = nn.Linear(2048, num_classes)
        for m in self.modules():
            if isinstance(m, nn.Conv2d) or isinstance(m, nn.
Linear):
                import scipy.stats as stats
                stddev = m.stddev if hasattr(m, 'stddev') else
0.1 #if else
                X = stats.truncnorm(-2, 2, scale=stddev)
                values = torch.Tensor(X.rvs(m.weight.data.nu
mel()))
                values = values.view(m.weight.data.size())
                m.weight.data.copy_(values)
            elif isinstance(m, nn.BatchNorm2d):
                m.weight.data.fill_(1)
                m.bias.data.zero_()

    def forward(self, x):
        x = x.view(-1, 3, INPUT_SIZE[0], INPUT_SIZE[1])
        x = self.Conv2d_1a_3x3(x)
        x = self.Conv2d_2a_3x3(x)
        x = self.Conv2d_2b_3x3(x)
        x = F.max_pool2d(x, kernel_size=3, stride=2)
        x = self.Conv2d_3b_1x1(x)
        x = self.Conv2d_4a_3x3(x)
        x = F.max_pool2d(x, kernel_size=3, stride=2)
        x = self.Mixed_5b(x)
        x = self.Mixed_5c(x)
        x = self.Mixed_5d(x)
        x = self.Mixed_6a(x)
        x = self.Mixed_6b(x)
        x = self.Mixed_6c(x)
        x = self.Mixed_6d(x)
        x = self.Mixed_6e(x)
        x = self.Mixed_7a(x)
        x = self.Mixed_7b(x)
        x = self.Mixed_7c(x)
        x = F.avg_pool2d(x, kernel_size=8)
        x = F.dropout(x, training=self.training)
        x = x.view(x.size(0), -1)
        x = self.fc(x)
        return x
```

Model Code




# OSS

## Databricks

**(-) 阿里云** [www.aliyun.com](http://www.aliyun.com)

# Evaluation






Confusion Matrix




	Score
F1-Score	0.811718
Precision	0.86934
Recall	0.819014
Accuracy	0.819014

Metrics



ROC Curve



PR Curve

**Thank for Listening!**